

Creation Date 05-Oct-2015

Revision Date 22-Dec-2024

Revision Number 1

**Section 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING****1.1. Product identifier**

Product Description: Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)  
Cat No. : T001514000

**1.2. Relevant identified uses of the substance or mixture and uses advised against**

Recommended Use Laboratory chemicals.  
Uses advised against No Information available

**1.3. Details of the supplier of the safety data sheet****Company**

**EU entity/business name**  
Thermo Fisher Scientific  
Janssen Pharmaceuticaaan 3a  
2440 Geel, Belgium

**UK entity/business name**  
Fisher Scientific UK  
Bishop Meadow Road, Loughborough,  
Leicestershire LE11 5RG, United Kingdom

**Swiss distributor** - Fisher Scientific AG  
Neuhofstrasse 11, CH 4153 Reinach  
Tel: +41 (0) 56 618 41 11  
e-mail - infoch@thermofisher.com

**E-mail address** begel.sdsdesk@thermofisher.com

**1.4. Emergency telephone number**

Tel: 01509 231166  
Chemtrec US: (800) 424-9300  
Chemtrec EU: 001-703-527-3887

For customers in Switzerland:  
Tox Info Suisse Emergency Number: **145 (24hr)**  
Tox Info Suisse: +41-44 251 51 51 (Emergency number from abroad)  
Chemtrec (24h) Toll-Free: 0800 564 402  
Chemtrec Local: +41-43 508 20 11 (Zurich)

**Section 2: HAZARDS IDENTIFICATION****2.1. Classification of the substance or mixture**

CLP Classification - Regulation (EC) No 1272/2008

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

## Physical hazards

Flammable liquids

Category 2 (H225)

## Health hazards

Serious Eye Damage/Eye Irritation  
Specific target organ toxicity - (single exposure)

Category 2 (H319)  
Category 3 (H336)

## Environmental hazards

Based on available data, the classification criteria are not met

Full text of Hazard Statements: see section 16

## 2.2. Label elements



Signal Word

Danger

## Hazard Statements

H225 - Highly flammable liquid and vapor  
H319 - Causes serious eye irritation  
H336 - May cause drowsiness or dizziness  
EUH066 - Repeated exposure may cause skin dryness or cracking

## Precautionary Statements

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray  
P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing  
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
P312 - Call a POISON CENTER or doctor if you feel unwell  
P403 + P233 - Store in a well-ventilated place. Keep container tightly closed  
P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking

## 2.3. Other hazards

This product does not contain any known or suspected endocrine disruptors

## Section 3: Composition/information on ingredients

### 3.2. Mixtures

| Component     | CAS No   | EC No             | Weight % | CLP Classification - Regulation (EC) No 1272/2008 |
|---------------|----------|-------------------|----------|---------------------------------------------------|
| Ethyl acetate | 141-78-6 | EEC No. 205-500-4 | 74-76    | Flam. Liq. 2 (H225)                               |

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

|               |         |           |       |                                                   |
|---------------|---------|-----------|-------|---------------------------------------------------|
|               |         |           |       | Eye Irrit. 2 (H319)<br>STOT SE 3 (H336)<br>EUH066 |
| Ethyl alcohol | 64-17-5 | 200-578-6 | 24-26 | Flam. Liq. 2 (H225)<br>Eye Irrit. 2 (H319)        |

| Component     | Specific concentration limits (SCL's) | M-Factor | Component notes |
|---------------|---------------------------------------|----------|-----------------|
| Ethyl alcohol | Eye Irrit. 2 :: C>=50%                | -        | -               |

Full text of Hazard Statements: see section 16

## Section 4: First aid measures

### 4.1. Description of first aid measures

|                                           |                                                                                                                                                  |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Eye Contact</b>                        | Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.                             |
| <b>Skin Contact</b>                       | Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes.                                            |
| <b>Ingestion</b>                          | Clean mouth with water and drink afterwards plenty of water.                                                                                     |
| <b>Inhalation</b>                         | Remove to fresh air.                                                                                                                             |
| <b>Self-Protection of the First Aider</b> | Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. |

### 4.2. Most important symptoms and effects, both acute and delayed

Difficulty in breathing. . Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting: Inhalation of high vapor concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting

### 4.3. Indication of any immediate medical attention and special treatment needed

|                           |                                                 |
|---------------------------|-------------------------------------------------|
| <b>Notes to Physician</b> | Treat symptomatically. Symptoms may be delayed. |
|---------------------------|-------------------------------------------------|

## Section 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable Extinguishing Media

Water mist may be used to cool closed containers.

#### Extinguishing media which must not be used for safety reasons

Do not use a solid water stream as it may scatter and spread fire.

### 5.2. Special hazards arising from the substance or mixture

Flammable. Containers may explode when heated. Vapors may form explosive mixtures with air. Vapors may travel to source of ignition and flash back.

#### Hazardous Combustion Products

None under normal use conditions.

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

## 5.3. Advice for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## Section 6: ACCIDENTAL RELEASE MEASURES

### 6.1. Personal precautions, protective equipment and emergency procedures

Remove all sources of ignition. Take precautionary measures against static discharges.

### 6.2. Environmental precautions

Do not flush into surface water or sanitary sewer system. See Section 12 for additional Ecological Information.

### 6.3. Methods and material for containment and cleaning up

Remove all sources of ignition. Use spark-proof tools and explosion-proof equipment. Take precautionary measures against static discharges.

### 6.4. Reference to other sections

Refer to protective measures listed in Sections 8 and 13.

## Section 7: Handling and storage

### 7.1. Precautions for safe handling

Keep away from open flames, hot surfaces and sources of ignition. Use only non-sparking tools. To avoid ignition of vapors by static electricity discharge, all metal parts of the equipment must be grounded. Take precautionary measures against static discharges. Use spark-proof tools and explosion-proof equipment.

#### **Hygiene Measures**

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Wash hands before breaks and after work.

### 7.2. Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place. Keep away from heat, sparks and flame. Flammables area.

**Technical Rules for Hazardous Substances (TRGS) 510**  
**Storage Class (LGK) (Germany)**

Class 3

**Switzerland - Storage of hazardous substances**

Storage class - SC 3  
<https://www.kvu.ch/de/themen/stoffe-und-produkte>  
<https://www.kvu.ch/fr/themes/substances-et-produits>  
<https://www.kvu.ch/it/temi/sostanze-e-prodotti>

### 7.3. Specific end use(s)

Use in laboratories

## Section 8: Exposure controls/personal protection

### 8.1. Control parameters

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

## Exposure limits

List source(s): **UK** - EH40/2005 Work Exposure Limits, Forth edition. Published 2020. **IRE** - 2010 Code of Practice for the Safety, Health and Welfare at Work (Chemical Agents) Regulations 2001. Published by the Health and Safety Authority. **CH** - The Government of Switzerland has set a directive on limit values for working materials (Grenzwerte am Arbeitsplatz) which is based on the Swiss Federal Regulation "Verordnung über die Verhütung von Unfällen und Berufskrankheiten". This directive is administered, periodically revised and enforced by SUVA (Swiss National Accident Insurance Fund). **EU** - Commission Directive (EU) 2019/1831 of 24 October 2019 establishing a fifth list of indicative occupational exposure limit values pursuant to Council Directive 98/24/EC and amending Commission Directive 2000/39/EC

| Component     | European Union                                                                                                        | The United Kingdom                                                                                                  | France                                                                                                                                                                                                                   | Belgium                                                                                                                         | Spain                                                                                                                                                                           |
|---------------|-----------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Ethyl acetate | TWA: 734 mg/m <sup>3</sup> (8h)<br>TWA: 200 ppm (8h)<br>STEL: 1468 mg/m <sup>3</sup> (15min)<br>STEL: 400 ppm (15min) | STEL: 1468 mg/m <sup>3</sup> 15 min<br>STEL: 400 ppm 15 min<br>TWA: 734 mg/m <sup>3</sup> 8 hr<br>TWA: 200 ppm 8 hr | TWA / VME: 200 ppm (8 heures). restrictive limit<br>TWA / VME: 734 mg/m <sup>3</sup> (8 heures). restrictive limit<br>STEL / VLCT: 400 ppm. restrictive limit<br>STEL / VLCT: 1468 mg/m <sup>3</sup> . restrictive limit | TWA: 200 ppm 8 uren<br>TWA: 734 mg/m <sup>3</sup> 8 uren<br>STEL: 400 ppm 15 minuten<br>STEL: 1468 mg/m <sup>3</sup> 15 minuten | STEL / VLA-EC: 400 ppm (15 minutos).<br>STEL / VLA-EC: 1468 mg/m <sup>3</sup> (15 minutos).<br>TWA / VLA-ED: 200 ppm (8 horas)<br>TWA / VLA-ED: 734 mg/m <sup>3</sup> (8 horas) |
| Ethyl alcohol |                                                                                                                       | TWA: 1000 ppm TWA;<br>1920 mg/m <sup>3</sup> TWA<br>WEL - STEL: 3000 ppm<br>STEL: 5760 mg/m <sup>3</sup> STEL       | TWA / VME: 1000 ppm (8 heures).<br>TWA / VME: 1900 mg/m <sup>3</sup> (8 heures).<br>STEL / VLCT: 5000 ppm.<br>STEL / VLCT: 9500 mg/m <sup>3</sup> .                                                                      | TWA: 1000 ppm 8 uren<br>TWA: 1907 mg/m <sup>3</sup> 8 uren                                                                      | STEL / VLA-EC: 1000 ppm (15 minutos).<br>STEL / VLA-EC: 1910 mg/m <sup>3</sup> (15 minutos).                                                                                    |

| Component     | Italy                                                                                                                                                                                                   | Germany                                                                                                                                                                                                                                                         | Portugal                                                                                                                          | The Netherlands                                                                                                                          | Finland                                                                                                                                            |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Ethyl acetate | TWA: 734 mg/m <sup>3</sup> 8 ore.<br>Time Weighted Average<br>TWA: 200 ppm 8 ore.<br>Time Weighted Average<br>STEL: 1468 mg/m <sup>3</sup> 15 minuti. Short-term<br>STEL: 400 ppm 15 minuti. Short-term | TWA: 200 ppm (8 Stunden). AGW - exposure factor 2<br>TWA: 730 mg/m <sup>3</sup> (8 Stunden). AGW - exposure factor 2<br>TWA: 200 ppm (8 Stunden). MAK<br>TWA: 750 mg/m <sup>3</sup> (8 Stunden). MAK<br>Höhepunkt: 400 ppm<br>Höhepunkt: 1500 mg/m <sup>3</sup> | STEL: 1468 mg/m <sup>3</sup> 15 minutos<br>STEL: 400 ppm 15 minutos<br>TWA: 200 ppm 8 horas<br>TWA: 734 mg/m <sup>3</sup> 8 horas | STEL: 400 ppm 15 minuten<br>STEL: 1468 mg/m <sup>3</sup> 15 minuten<br>TWA: 200 ppm 8 uren<br>TWA: 734 mg/m <sup>3</sup> 8 uren          | TWA: 200 ppm 8 tunteina<br>TWA: 730 mg/m <sup>3</sup> 8 tunteina<br>STEL: 400 ppm 15 minuutteina<br>STEL: 1470 mg/m <sup>3</sup> 15 minuutteina    |
| Ethyl alcohol |                                                                                                                                                                                                         | 200 ppm TWA MAK;<br>380 mg/m <sup>3</sup> TWA MAK                                                                                                                                                                                                               | STEL: 1000 ppm 15 minutos                                                                                                         | huid<br>STEL: 1000 ppm 15 minuten<br>STEL: 1900 mg/m <sup>3</sup> 15 minuten<br>TWA: 137 ppm 8 uren<br>TWA: 260 mg/m <sup>3</sup> 8 uren | TWA: 1000 ppm 8 tunteina<br>TWA: 1900 mg/m <sup>3</sup> 8 tunteina<br>STEL: 1300 ppm 15 minuutteina<br>STEL: 2500 mg/m <sup>3</sup> 15 minuutteina |

| Component     | Austria                                                                                                                                                  | Denmark                                                                                                                                | Switzerland                                                                                                                            | Poland                                                                             | Norway                                                                                                                                                                                    |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Ethyl acetate | MAK-KZGW: 400 ppm 15 Minuten<br>MAK-KZGW: 1468 mg/m <sup>3</sup> 15 Minuten<br>MAK-TMW: 200 ppm 8 Stunden<br>MAK-TMW: 734 mg/m <sup>3</sup> 8 Stunden    | TWA: 150 ppm 8 timer<br>TWA: 540 mg/m <sup>3</sup> 8 timer<br>STEL: 1468 mg/m <sup>3</sup> 15 minutter<br>STEL: 400 ppm 15 minutter    | STEL: 400 ppm 15 Minuten<br>STEL: 1460 mg/m <sup>3</sup> 15 Minuten<br>TWA: 200 ppm 8 Stunden<br>TWA: 730 mg/m <sup>3</sup> 8 Stunden  | STEL: 1468 mg/m <sup>3</sup> 15 minutach<br>TWA: 734 mg/m <sup>3</sup> 8 godzinach | TWA: 200 ppm 8 timer<br>TWA: 734 mg/m <sup>3</sup> 8 timer<br>STEL: 400 ppm 15 minutter. value from the regulation<br>STEL: 1468 mg/m <sup>3</sup> 15 minutter. value from the regulation |
| Ethyl alcohol | MAK-KZGW: 2000 ppm 15 Minuten<br>MAK-KZGW: 3800 mg/m <sup>3</sup> 15 Minuten<br>MAK-TMW: 1000 ppm 8 Stunden<br>MAK-TMW: 1900 mg/m <sup>3</sup> 8 Stunden | TWA: 1000 ppm 8 timer<br>TWA: 1900 mg/m <sup>3</sup> 8 timer<br>STEL: 2000 ppm 15 minutter<br>STEL: 3800 mg/m <sup>3</sup> 15 minutter | STEL: 1000 ppm 15 Minuten<br>STEL: 1920 mg/m <sup>3</sup> 15 Minuten<br>TWA: 500 ppm 8 Stunden<br>TWA: 960 mg/m <sup>3</sup> 8 Stunden | TWA: 1900 mg/m <sup>3</sup> 8 godzinach                                            | TWA: 500 ppm 8 timer<br>TWA: 950 mg/m <sup>3</sup> 8 timer<br>STEL: 625 ppm 15 minutter. value calculated<br>STEL: 1187.5 mg/m <sup>3</sup> 15 minutter. value calculated                 |

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

| Component     | Bulgaria                                                                                      | Croatia                                                                                                                                                     | Ireland                                                                                                               | Cyprus                                                                                      | Czech Republic                                                             |
|---------------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|
| Ethyl acetate | TWA: 734 mg/m <sup>3</sup><br>TWA: 200 ppm<br>STEL : 1468 mg/m <sup>3</sup><br>STEL : 400 ppm | TWA-GVI: 200 ppm 8 satima.<br>TWA-GVI: 734 mg/m <sup>3</sup> 8 satima.<br>STEL-KGVI: 400 ppm 15 minutama.<br>STEL-KGVI: 1468 mg/m <sup>3</sup> 15 minutama. | TWA: 734 mg/m <sup>3</sup> 8 hr.<br>TWA: 200 ppm 8 hr.<br>STEL: 1468 mg/m <sup>3</sup> 15 min<br>STEL: 400 ppm 15 min | STEL: 1468 mg/m <sup>3</sup><br>STEL: 400 ppm<br>TWA: 734 mg/m <sup>3</sup><br>TWA: 200 ppm | TWA: 700 mg/m <sup>3</sup> 8 hodinách.<br>Ceiling: 900 mg/m <sup>3</sup>   |
| Ethyl alcohol | TWA: 1000 mg/m <sup>3</sup>                                                                   | TWA-GVI: 1000 ppm 8 satima.<br>TWA-GVI: 1900 mg/m <sup>3</sup> 8 satima.                                                                                    | STEL: 1000 ppm 15 min                                                                                                 |                                                                                             | TWA: 1000 mg/m <sup>3</sup> 8 hodinách.<br>Ceiling: 3000 mg/m <sup>3</sup> |

| Component     | Estonia                                                                                                                                           | Gibraltar                                                                                                           | Greece                                                                                      | Hungary                                                                                                                                                    | Iceland                                                                                                                                |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Ethyl acetate | TWA: 150 ppm 8 tundides.<br>TWA: 500 mg/m <sup>3</sup> 8 tundides.<br>STEL: 300 ppm 15 minutites.<br>STEL: 1100 mg/m <sup>3</sup> 15 minutites.   | TWA: 734 ppm 8 hr<br>TWA: 200 mg/m <sup>3</sup> 8 hr<br>STEL: 1468 ppm 15 min<br>STEL: 400 mg/m <sup>3</sup> 15 min | STEL: 400 ppm<br>STEL: 1468 mg/m <sup>3</sup><br>TWA: 200 ppm<br>TWA: 734 mg/m <sup>3</sup> | STEL: 400 ppm 15 percekben. CK<br>STEL: 1468 mg/m <sup>3</sup> 15 percekben. CK<br>TWA: 200 ppm 8 órában. AK<br>TWA: 734 mg/m <sup>3</sup> 8 órában. AK    | TWA: 150 ppm 8 klukkustundum.<br>TWA: 540 mg/m <sup>3</sup> 8 klukkustundum.<br>Ceiling: 300 ppm<br>Ceiling: 1080 mg/m <sup>3</sup>    |
| Ethyl alcohol | TWA: 500 ppm 8 tundides.<br>TWA: 1000 mg/m <sup>3</sup> 8 tundides.<br>STEL: 1000 ppm 15 minutites.<br>STEL: 1900 mg/m <sup>3</sup> 15 minutites. |                                                                                                                     | TWA: 1000 ppm<br>TWA: 1900 mg/m <sup>3</sup>                                                | STEL: 2000 ppm 15 percekben. CK<br>STEL: 3800 mg/m <sup>3</sup> 15 percekben. CK<br>TWA: 1000 ppm 8 órában. AK<br>TWA: 1900 mg/m <sup>3</sup> 8 órában. AK | TWA: 1000 ppm 8 klukkustundum.<br>TWA: 1900 mg/m <sup>3</sup> 8 klukkustundum.<br>Ceiling: 2000 ppm<br>Ceiling: 3800 mg/m <sup>3</sup> |

| Component     | Latvia                                                                                     | Lithuania                                                                                                   | Luxembourg                                                                                                                            | Malta                                                                                                           | Romania                                                                                                                        |
|---------------|--------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Ethyl acetate | STEL: 1468 mg/m <sup>3</sup><br>STEL: 400 ppm<br>TWA: 200 mg/m <sup>3</sup><br>TWA: 54 ppm | Ceiling: 300 ppm<br>Ceiling: 1100 mg/m <sup>3</sup><br>TWA: 150 ppm IPRD<br>TWA: 500 mg/m <sup>3</sup> IPRD | TWA: 734 mg/m <sup>3</sup> 8 Stunden<br>TWA: 200 ppm 8 Stunden<br>STEL: 1468 mg/m <sup>3</sup> 15 Minuten<br>STEL: 400 ppm 15 Minuten | TWA: 200 ppm<br>TWA: 734 mg/m <sup>3</sup><br>STEL: 400 ppm 15 minuti<br>STEL: 1468 mg/m <sup>3</sup> 15 minuti | TWA: 200 ppm 8 ore<br>TWA: 734 mg/m <sup>3</sup> 8 ore<br>STEL: 400 ppm 15 minute<br>STEL: 1468 mg/m <sup>3</sup> 15 minute    |
| Ethyl alcohol | TWA: 1000 mg/m <sup>3</sup>                                                                | TWA: 500 ppm IPRD<br>TWA: 1000 mg/m <sup>3</sup> IPRD<br>STEL: 1000 ppm<br>STEL: 1900 mg/m <sup>3</sup>     |                                                                                                                                       |                                                                                                                 | TWA: 1000 ppm 8 ore<br>TWA: 1900 mg/m <sup>3</sup> 8 ore<br>STEL: 5000 ppm 15 minute<br>STEL: 9500 mg/m <sup>3</sup> 15 minute |

| Component     | Russia                                                          | Slovak Republic                                                               | Slovenia                                                                                                                         | Sweden                                                                                                                                                                | Turkey |
|---------------|-----------------------------------------------------------------|-------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| Ethyl acetate | TWA: 50 mg/m <sup>3</sup> 2417<br>MAC: 200 mg/m <sup>3</sup>    | Ceiling: 1100 mg/m <sup>3</sup><br>TWA: 200 ppm<br>TWA: 734 mg/m <sup>3</sup> | TWA: 200 ppm 8 urah<br>TWA: 734 mg/m <sup>3</sup> 8 urah<br>STEL: 400 ppm 15 minutah<br>STEL: 1468 mg/m <sup>3</sup> 15 minutah  | Binding STEL: 300 ppm 15 minuter<br>Binding STEL: 1100 mg/m <sup>3</sup> 15 minuter<br>TLV: 150 ppm 8 timmar. NGV<br>TLV: 550 mg/m <sup>3</sup> 8 timmar. NGV         |        |
| Ethyl alcohol | TWA: 1000 mg/m <sup>3</sup> 2391<br>MAC: 2000 mg/m <sup>3</sup> | Ceiling: 1920 mg/m <sup>3</sup><br>TWA: 500 ppm<br>TWA: 960 mg/m <sup>3</sup> | TWA: 960 mg/m <sup>3</sup> 8 urah<br>TWA: 500 ppm 8 urah<br>STEL: 1000 ppm 15 minutah<br>STEL: 1920 mg/m <sup>3</sup> 15 minutah | Indicative STEL: 1000 ppm 15 minuter<br>Indicative STEL: 1900 mg/m <sup>3</sup> 15 minuter<br>TLV: 500 ppm 8 timmar. NGV<br>TLV: 1000 mg/m <sup>3</sup> 8 timmar. NGV |        |

## Biological limit values

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

## Monitoring methods

BS EN 14042:2003 Title Identifier: Workplace atmospheres. Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents.

MDHS70 General methods for sampling airborne gases and vapours

MDHS 88 Volatile organic compounds in air. Laboratory method using diffusive samplers, solvent desorption and gas chromatography

MDHS 96 Volatile organic compounds in air - Laboratory method using pumped solid sorbent tubes, solvent desorption and gas chromatography

## Derived No Effect Level (DNEL) / Derived Minimum Effect Level (DMEL)

See table for values

| Component                          | Acute effects local (Oral) | Acute effects systemic (Oral) | Chronic effects local (Oral) | Chronic effects systemic (Oral) |
|------------------------------------|----------------------------|-------------------------------|------------------------------|---------------------------------|
| Ethyl alcohol<br>64-17-5 ( 24-26 ) |                            | DNEL = 87 mg/kg bw/d          |                              |                                 |

| Component                           | Acute effects local (Dermal) | Acute effects systemic (Dermal) | Chronic effects local (Dermal) | Chronic effects systemic (Dermal) |
|-------------------------------------|------------------------------|---------------------------------|--------------------------------|-----------------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) |                              |                                 |                                | DNEL = 63mg/kg bw/day             |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  |                              |                                 |                                | DNEL = 343mg/kg bw/day            |

| Component                           | Acute effects local (Inhalation)         | Acute effects systemic (Inhalation)      | Chronic effects local (Inhalation)      | Chronic effects systemic (Inhalation) |
|-------------------------------------|------------------------------------------|------------------------------------------|-----------------------------------------|---------------------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | DNEL = 1468 mg/m <sup>3</sup><br>400 ppm | DNEL = 1468 mg/m <sup>3</sup><br>400 ppm | DNEL = 734 mg/m <sup>3</sup><br>200 ppm | DNEL = 734mg/m <sup>3</sup>           |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  | DNEL = 1900mg/m <sup>3</sup>             |                                          |                                         | DNEL = 950mg/m <sup>3</sup>           |

## Predicted No Effect Concentration (PNEC)

See values below.

| Component                           | Fresh water     | Fresh water sediment         | Water Intermittent | Microorganisms in sewage treatment | Soil (Agriculture)        |
|-------------------------------------|-----------------|------------------------------|--------------------|------------------------------------|---------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | PNEC = 0.24mg/L | PNEC = 1.15mg/kg sediment dw | PNEC = 1.65mg/L    | PNEC = 650mg/L                     | PNEC = 0.148mg/kg soil dw |

| Component                           | Marine water     | Marine water sediment         | Marine water Intermittent | Food chain          | Air |
|-------------------------------------|------------------|-------------------------------|---------------------------|---------------------|-----|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | PNEC = 0.024mg/L | PNEC = 0.115mg/kg sediment dw |                           | PNEC = 0.2g/kg food |     |

## 8.2. Exposure controls

### Engineering Measures

Ensure that eyewash stations and safety showers are close to the workstation location. Ensure adequate ventilation, especially in confined areas. Use explosion-proof electrical/ventilating/lighting/equipment.

Wherever possible, engineering control measures such as the isolation or enclosure of the process, the introduction of process or

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

equipment changes to minimise release or contact, and the use of properly designed ventilation systems, should be adopted to control hazardous materials at source

## Personal protective equipment

**Eye Protection** Goggles (European standard - EN 166)

**Hand Protection** Protective gloves

| Glove material              | Breakthrough time                    | Glove thickness | EU standard | Glove comments        |
|-----------------------------|--------------------------------------|-----------------|-------------|-----------------------|
| Nitrile rubber<br>Viton (R) | See manufacturers<br>recommendations | -               | EN 374      | (minimum requirement) |

**Skin and body protection** Long sleeved clothing.

Inspect gloves before use. observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information) gloves are suitable for the task: Chemical compatability, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion. gloves with care avoiding skin contamination.

**Respiratory Protection** When workers are facing concentrations above the exposure limit they must use appropriate certified respirators.  
To protect the wearer, respiratory protective equipment must be the correct fit and be used and maintained properly

**Large scale/emergency use** Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced  
**Recommended Filter type:** low boiling organic solvent Type AX Brown conforming to EN371 or Organic gases and vapours filter Type A Brown conforming to EN14387

**Small scale/Laboratory use** Use a NIOSH/MSHA or European Standard EN 149:2001 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.  
**Recommended half mask:-** Valve filtering: EN405; or; Half mask: EN140; plus filter, EN 141  
When RPE is used a face piece Fit Test should be conducted

**Environmental exposure controls** Prevent product from entering drains. Do not allow material to contaminate ground water system.

## Section 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

|                                                |                          |                                               |
|------------------------------------------------|--------------------------|-----------------------------------------------|
| <b>Physical State</b>                          | Liquid                   |                                               |
| <b>Appearance</b>                              | Colorless                |                                               |
| <b>Odor</b>                                    | No information available |                                               |
| <b>Odor Threshold</b>                          | No data available        |                                               |
| <b>Melting Point/Range</b>                     | No data available        |                                               |
| <b>Softening Point</b>                         | No data available        |                                               |
| <b>Boiling Point/Range</b>                     | No information available |                                               |
| <b>Flammability (liquid)</b>                   | Highly flammable         | On basis of test data                         |
| <b>Flammability (solid,gas)</b>                | Not applicable           | Liquid                                        |
| <b>Explosion Limits</b>                        | No data available        |                                               |
| <b>Flash Point</b>                             | -5 °C / 23 °F            | <b>Method -</b> Based on available literature |
| <b>Autoignition Temperature</b>                | No data available        |                                               |
| <b>Decomposition Temperature</b>               | No data available        |                                               |
| <b>pH</b>                                      | No information available |                                               |
| <b>Viscosity</b>                               | No data available        |                                               |
| <b>Water Solubility</b>                        | Miscible                 |                                               |
| <b>Solubility in other solvents</b>            | No information available |                                               |
| <b>Partition Coefficient (n-octanol/water)</b> |                          |                                               |

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

|                                   |                         |             |
|-----------------------------------|-------------------------|-------------|
| <b>Component</b>                  | <b>log Pow</b>          |             |
| Ethyl acetate                     | 0.73                    |             |
| Ethyl alcohol                     | -0.32                   |             |
| <b>Vapor Pressure</b>             | No data available       |             |
| <b>Density / Specific Gravity</b> | 0.87                    | Estimated   |
| <b>Bulk Density</b>               | Not applicable          | Liquid      |
| <b>Vapor Density</b>              | No data available       | (Air = 1.0) |
| <b>Particle characteristics</b>   | Not applicable (liquid) |             |

## 9.2. Other information

**Explosive Properties** Vapors may form explosive mixtures with air

## Section 10: Stability and reactivity

**10.1. Reactivity** None known, based on information available

**10.2. Chemical stability** Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

**Hazardous Polymerization** No information available.  
**Hazardous Reactions** No information available.

**10.4. Conditions to avoid** Keep away from open flames, hot surfaces and sources of ignition.

**10.5. Incompatible materials** Strong oxidizing agents.

**10.6. Hazardous decomposition products** None under normal use conditions.

## Section 11: Toxicological information

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Product Information

**(a) acute toxicity;**  
**Oral** Based on available data, the classification criteria are not met  
**Dermal** Based on available data, the classification criteria are not met  
**Inhalation** Based on available data, the classification criteria are not met

#### Toxicology data for the components

| Component     | LD50 Oral                                                    | LD50 Dermal                                       | LC50 Inhalation                                                   |
|---------------|--------------------------------------------------------------|---------------------------------------------------|-------------------------------------------------------------------|
| Ethyl acetate | 10,200 mg/kg ( Rat )                                         | > 20 mL/kg ( Rabbit )<br>> 18000 mg/kg ( Rabbit ) | 58 mg/l (rat; 8 h)                                                |
| Ethyl alcohol | LD50 = 10470 mg/kg<br>OECD 401 (Rat)<br>3450 mg/kg ( Mouse ) | -                                                 | LC50 = 117-125 mg/l (4h)<br>OECD 403 (rat)<br>20000 ppm/10H (rat) |

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

(b) skin corrosion/irritation; No data available

(c) serious eye damage/irritation; Category 2

(d) respiratory or skin sensitization;

Respiratory No data available  
Skin No data available

| Component                           | Test method                                       | Test species | Study result      |
|-------------------------------------|---------------------------------------------------|--------------|-------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | OECD Test Guideline 406                           | guinea pig   | - non-sensitising |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  | Mouse Ear Swelling Test (MEST)                    | mouse        | non-sensitising   |
|                                     | OECD Test Guideline 429<br>Local Lymph Node Assay | mouse        | non-sensitising   |

(e) germ cell mutagenicity; No data available

| Component                           | Test method                                             | Test species          | Study result |
|-------------------------------------|---------------------------------------------------------|-----------------------|--------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | OECD Test Guideline 471<br>AMES test                    | in vitro<br>Bacteria  | negative     |
|                                     | OECD Test Guideline 473<br>Chromosomal aberration assay | in vitro<br>Mammalian | negative     |
|                                     | OECD Test Guideline 476<br>Gene cell mutation           | in vitro<br>Mammalian | negative     |
|                                     | OECD Test Guideline 474<br>Mouse micronucleus assay     | in vivo<br>Mammalian  | negative     |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  | AMES test<br>OECD Test Guideline 471                    | in vitro<br>Bacteria  | negative     |
|                                     | Gene cell mutation<br>OECD Test Guideline 476           | in vitro<br>Mammalian | negative     |

(f) carcinogenicity; No data available

The table below indicates whether each agency has listed any ingredient as a carcinogen

(g) reproductive toxicity; No data available

| Component                           | Test method             | Test species / Duration       | Study result                       |
|-------------------------------------|-------------------------|-------------------------------|------------------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | OECD Test Guideline 416 | Oral<br>mouse<br>2 Generation | NOAEL =<br>26400<br>mg/kg bw/day   |
|                                     | OECD Test Guideline 414 | Inhalation<br>Rat             | NOAEC =<br>73300 mg/m <sup>3</sup> |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  | OECD Test Guideline 416 | Oral / mouse<br>2 Generation  | NOAEL = 13.8 g/kg/day              |
|                                     | OECD Test Guideline 414 | Inhalation / Rat              | NOAEC =<br>16000 ppm               |

(h) STOT-single exposure; Category 3

Results / Target organs Central nervous system (CNS).

(i) STOT-repeated exposure; No data available

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

**Target Organs** No information available.

**(j) aspiration hazard;** No data available

**Symptoms / effects, both acute and delayed** Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting. Inhalation of high vapor concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.

## 11.2. Information on other hazards

**Endocrine Disrupting Properties** Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

## Section 12: Ecological information

### 12.1. Toxicity

**Ecotoxicity effects** Contains no substances known to be hazardous to the environment or that are not degradable in waste water treatment plants.

| Component     | Freshwater Fish                                                      | Water Flea                                    | Freshwater Algae                           |
|---------------|----------------------------------------------------------------------|-----------------------------------------------|--------------------------------------------|
| Ethyl acetate | Fathead minnow: LC50: 230 mg/l/ 96h<br>Gold orfe: LC50: 270 mg/L/48h | EC50 = 717 mg/L/48h                           | EC50 = 3300 mg/L/48h                       |
| Ethyl alcohol | Fathead minnow (Pimephales promelas) LC50 = 14200 mg/l/96h           | EC50 = 9268 mg/L/48h<br>EC50 = 10800 mg/L/24h | EC50 (72h) = 275 mg/l (Chlorella vulgaris) |

| Component     | Microtox                                                                                                    | M-Factor |
|---------------|-------------------------------------------------------------------------------------------------------------|----------|
| Ethyl acetate | EC50 = 1180 mg/L 5 min<br>EC50 = 1500 mg/L 15 min<br>EC50 = 5870 mg/L 15 min<br>EC50 = 7400 mg/L 2 h        |          |
| Ethyl alcohol | Photobacterium phosphoreum: EC50 = 34634 mg/L/30 min<br>Photobacterium phosphoreum: EC50 = 35470 mg/L/5 min |          |

### 12.2. Persistence and degradability

**Persistence** Miscible with water, Persistence is unlikely, based on information available.

| Component                           | Degradability            |
|-------------------------------------|--------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) | 79 % (20 d) (OECD 301 D) |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  | OECD 301E = 94%          |

**Degradation in sewage treatment plant** Contains substances known to be hazardous to the environment or not degradable in waste water treatment plants.

**12.3. Bioaccumulative potential** Bioaccumulation is unlikely

| Component     | log Pow | Bioconcentration factor (BCF) |
|---------------|---------|-------------------------------|
| Ethyl acetate | 0.73    | 30 dimensionless              |
| Ethyl alcohol | -0.32   | No data available             |

**12.4. Mobility in soil** The product is water soluble, and may spread in water systems Will likely be mobile in the environment due to its volatility. Highly mobile in soils

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

**12.5. Results of PBT and vPvB assessment** No data available for assessment.

**12.6. Endocrine disrupting properties**

**Endocrine Disruptor Information** This product does not contain any known or suspected endocrine disruptors

**12.7. Other adverse effects**  
**Persistent Organic Pollutant**  
**Ozone Depletion Potential**

This product does not contain any known or suspected substance  
This product does not contain any known or suspected substance

## Section 13: Disposal considerations

### 13.1. Waste treatment methods

**Waste from Residues/Unused Products**

Waste is classified as hazardous. Dispose of in accordance with the European Directives on waste and hazardous waste. Dispose of in accordance with local regulations.

**Contaminated Packaging**

Dispose of this container to hazardous or special waste collection point. Empty containers retain product residue, (liquid and/or vapor), and can be dangerous. Keep product and empty container away from heat and sources of ignition.

**European Waste Catalogue (EWC)**

According to the European Waste Catalog, Waste Codes are not product specific, but application specific.

**Other Information**

Do not flush to sewer. Waste codes should be assigned by the user based on the application for which the product was used. Can be landfilled or incinerated, when in compliance with local regulations.

**Switzerland - Waste Ordinance**

Disposal should be in accordance with applicable regional, national and local laws and regulations. Ordinance on the Avoidance and the Disposal of Waste (Waste Ordinance, ADWO) SR 814.600  
<https://www.fedlex.admin.ch/eli/cc/2015/891/en>

## Section 14: Transport information

### IMDG/IMO

**14.1. UN number**

UN1993

**14.2. UN proper shipping name**  
**Technical Shipping Name**

Flammable liquid, n.o.s.  
Ethylacetate, Ethyl alcohol

**14.3. Transport hazard class(es)**

3

**14.4. Packing group**

II

### ADR

**14.1. UN number**

UN1993

**14.2. UN proper shipping name**  
**Technical Shipping Name**

Flammable liquid, n.o.s.  
Ethylacetate, Ethyl alcohol

**14.3. Transport hazard class(es)**

3

**14.4. Packing group**

II

### IATA

**14.1. UN number**

UN1993

**14.2. UN proper shipping name**  
**Technical Shipping Name**

Flammable liquid, n.o.s.  
Ethylacetate, Ethyl alcohol

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

**14.3. Transport hazard class(es)** 3

**14.4. Packing group** II

**14.5. Environmental hazards** No hazards identified

**14.6. Special precautions for user** No special precautions required.

**14.7. Maritime transport in bulk according to IMO instruments** Not applicable, packaged goods

## Section 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### International Inventories

Europe (EINECS/ELINCS/NLP), China (IECSC), Taiwan (TCSI), Korea (KECL), Japan (ENCS), Japan (ISHL), Canada (DSL/NDSL), Australia (AICS), New Zealand (NZIoC), Philippines (PICCS). US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

| Component     | CAS No   | EINECS    | ELINCS | NLP | IECSC | TCSI | KECL     | ENCS | ISHL |
|---------------|----------|-----------|--------|-----|-------|------|----------|------|------|
| Ethyl acetate | 141-78-6 | 205-500-4 | -      | -   | X     | X    | KE-00047 | X    | X    |
| Ethyl alcohol | 64-17-5  | 200-578-6 | -      | -   | X     | X    | KE-13217 | X    | X    |

| Component     | CAS No   | TSCA | TSCA Inventory notification - Active-Inactive | DSL | NDSL | AICS | NZIoC | PICCS |
|---------------|----------|------|-----------------------------------------------|-----|------|------|-------|-------|
| Ethyl acetate | 141-78-6 | X    | ACTIVE                                        | X   | -    | X    | X     | X     |
| Ethyl alcohol | 64-17-5  | X    | ACTIVE                                        | X   | -    | X    | X     | X     |

**Legend:** X - Listed '-' - Not Listed

**KECL** - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

#### Authorisation/Restrictions according to EU REACH

| Component     | CAS No   | REACH (1907/2006) - Annex XIV - Substances Subject to Authorization | REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances | REACH Regulation (EC 1907/2006) article 59 - Candidate List of Substances of Very High Concern (SVHC) |
|---------------|----------|---------------------------------------------------------------------|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| Ethyl acetate | 141-78-6 | -                                                                   | Use restricted. See entry 75. (see link for restriction details)              | -                                                                                                     |
| Ethyl alcohol | 64-17-5  | -                                                                   | -                                                                             | -                                                                                                     |

#### REACH links

<https://echa.europa.eu/substances-restricted-under-reach>

#### Seveso III Directive (2012/18/EC)

| Component     | CAS No   | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements |
|---------------|----------|-------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| Ethyl acetate | 141-78-6 | Not applicable                                                                            | Not applicable                                                                           |
| Ethyl alcohol | 64-17-5  | Not applicable                                                                            | Not applicable                                                                           |

Regulation (EC) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of dangerous chemicals

Not applicable

Contains component(s) that meet a 'definition' of per & poly fluoroalkyl substance (PFAS)?

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

Not applicable

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work .

Take note of Directive 2000/39/EC establishing a first list of indicative occupational exposure limit values

## National Regulations

**UK** - Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment

**WGK Classification** Water endangering class = 1 (self classification)

| Component     | Germany - Water Classification (AwSV) | Germany - TA-Luft Class |
|---------------|---------------------------------------|-------------------------|
| Ethyl acetate | WGK1                                  |                         |
| Ethyl alcohol | WGK1                                  |                         |

| Component     | France - INRS (Tables of occupational diseases)      |
|---------------|------------------------------------------------------|
| Ethyl acetate | Tableaux des maladies professionnelles (TMP) - RG 84 |
| Ethyl alcohol | Tableaux des maladies professionnelles (TMP) - RG 84 |

## Swiss Regulations

Article 4 para. 4 of the Ordinance on the protection of young people in the workplace (SR 822.115) and Article 1 lit. f of the EAER regulation on hazardous work and young people (SR 822.115.2).

Take note on Article 13 Maternity Ordinance (SR 822.111.52) with regards expectant and nursing mothers.

| Component                           | Switzerland - Ordinance on the Reduction of Risk from handling of hazardous substances preparation (SR 814.81) | Switzerland - Ordinance on Incentive Taxes on Volatile Organic Compounds (OVOC) | Switzerland - Ordinance of the Rotterdam Convention on the Prior Informed Consent Procedure |
|-------------------------------------|----------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| Ethyl acetate<br>141-78-6 ( 74-76 ) |                                                                                                                | Group I                                                                         |                                                                                             |
| Ethyl alcohol<br>64-17-5 ( 24-26 )  |                                                                                                                | Group I                                                                         |                                                                                             |

## 15.2. Chemical safety assessment

Chemical Safety Assessment/Reports (CSA/CSR) are not required for mixtures

## Section 16: Other information

### Full text of H-Statements referred to under sections 2 and 3

H225 - Highly flammable liquid and vapor

H319 - Causes serious eye irritation

H336 - May cause drowsiness or dizziness

EUH066 - Repeated exposure may cause skin dryness or cracking

### Legend

**CAS** - Chemical Abstracts Service

**EINECS/ELINCS** - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**IECSC** - Chinese Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**ENCS** - Japanese Existing and New Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**NZIoC** - New Zealand Inventory of Chemicals

# SAFETY DATA SHEET

Ethyl acetate / Ethanol Solution, for HPLC, 3:1 (v/v)

Revision Date 22-Dec-2024

**WEL** - Workplace Exposure Limit

**ACGIH** - American Conference of Governmental Industrial Hygienists

**DNEL** - Derived No Effect Level

**RPE** - Respiratory Protective Equipment

**LC50** - Lethal Concentration 50%

**NOEC** - No Observed Effect Concentration

**PBT** - Persistent, Bioaccumulative, Toxic

**TWA** - Time Weighted Average

**IARC** - International Agency for Research on Cancer

Predicted No Effect Concentration (PNEC)

**LD50** - Lethal Dose 50%

**EC50** - Effective Concentration 50%

**POW** - Partition coefficient Octanol:Water

**vPvB** - very Persistent, very Bioaccumulative

**ADR** - European Agreement Concerning the International Carriage of Dangerous Goods by Road

**IMO/IMDG** - International Maritime Organization/International Maritime Dangerous Goods Code

**OECD** - Organisation for Economic Co-operation and Development

**BCF** - Bioconcentration factor

**ICAO/IATA** - International Civil Aviation Organization/International Air Transport Association

**MARPOL** - International Convention for the Prevention of Pollution from Ships

**ATE** - Acute Toxicity Estimate

**VOC** - (volatile organic compound)

## Key literature references and sources for data

<https://echa.europa.eu/information-on-chemicals>

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

## Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

**Physical hazards** On basis of test data

**Health Hazards** Calculation method

**Environmental hazards** Calculation method

## Training Advice

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

Use of personal protective equipment, covering appropriate selection, compatibility, breakthrough thresholds, care, maintenance, fit and standards.

First aid for chemical exposure, including the use of eye wash and safety showers.

Fire prevention and fighting, identifying hazards and risks, static electricity, explosive atmospheres posed by vapours and dusts.

Chemical incident response training.

**Creation Date** 05-Oct-2015

**Revision Date** 22-Dec-2024

**Revision Summary** Not applicable.

**This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.  
COMMISSION REGULATION (EU) 2020/878 amending Annex II to Regulation (EC) No  
1907/2006 .**

**For Switzerland - Compiled in accordance with the technical provisions referred to in Annex 2,  
Number 3, ChemO (SR 813.11 - Ordinance on Protection against Dangerous Substances and  
Preparations).**

## Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet**