

Creation Date 07-Mar-2022

Revision Date 07-Mar-2022

Revision Number 1

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**1.1. Product identifier**

Product Code/Catalogue Number: 10028521, 10028522, 10028523, 10028524
Product Description: Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1, Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

1.2. Relevant identified uses of the substance or mixture and uses advised against

Recommended Use For Investigational Use/Performance Evaluation Only.
Uses advised against No Information available

1.3. Details of the supplier of the safety data sheet

Company Thermo Fisher Scientific Oy
Ratastie 2
FI-01620 Vantaa, Finland
+358 10 329200.
E-mail address cascadion.sds@thermofisher.com
www.thermofisher.com/cascadion
edfu.thermofisher.com

1.4. Emergency telephone number

CHEMTREC: +1-703-527-3887 (INTERNATIONAL)
1-800-424-9300 (NORTH AMERICA)

SECTION 2: HAZARDS IDENTIFICATION**2.1. Classification of the substance or mixture****CLP Classification - Regulation (EC) No 1272/2008****Physical hazards**

Based on available data, the classification criteria are not met

Health hazards

Based on available data, the classification criteria are not met

Environmental hazards

Based on available data, the classification criteria are not met

Full text of Hazard Statements: see section 16

2.2. Label elements

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

None required

EUH210 - Safety data sheet available on request

2.3. Other hazards

No information available

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.2. Mixtures

Component	CAS No	EC No	Weight %	CLP Classification - Regulation (EC) No 1272/2008
Human Serum	NA		>98	-
Ethyl alcohol	64-17-5	200-578-6	<1.2	Flam. Liq. 2 (H225)
Dimethyl sulfoxide	67-68-5	EEC No. 200-664-3	<0.6	-
Methyl alcohol	67-56-1	200-659-6	<0.2	Flam. Liq. 2 (H225) Acute Tox. 3 (H301) Acute Tox. 3 (H311) Acute Tox. 3 (H331) STOT SE 1 (H370)
Fluconazole	86386-73-4		<0.01	Acute Tox. 4 (H302) Skin Irrit. 2 (H315) Eye Irrit. 2 (H319) Repr. 1B (H360D)
2-Methyl-2H-isothiazol-3-one	2682-20-4	EEC No. 220-239-6	<0.01	Skin Corr. 1B (H314) Skin Sens. 1 (H317) Acute Tox. 4 (H302) Acute Tox. 3 (H331) Aquatic Acute 1 (H400)
Voriconazole	137234-62-9		<0.001	Acute Tox. 3 (H301) Carc. 2 (H351) Repr. 1A (H360) STOT RE 2 (H373) Aquatic Chronic 3 (H412)
Posaconazole	171228-49-2		<0.001	Repr. 2 (H361)
Isavuconazole	241479-67-4		<0.001	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-2-(2-hydroxy-1-methylpropyl)-	112559-91-8		<0.001	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-2-(1-methylpropyl)-	84625-61-6		<0.001	Acute Tox. 4 (H302) Eye Irrit. 2 (H319) Skin Irrit. 2 (H315) STOT SE 3 (H335)

Component	Specific concentration limits (SCL's)	M-Factor	Component notes
Methyl alcohol	STOT SE 1 (H370) :: C>=10% STOT SE 2 (H371) ::	-	-

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

	3%≤C<10%		
2-Methyl-2H-isothiazol-3-one	Skin Sens. 1A (H317) :: C>=0.0015%	10 (acute) 1 (chronic)	-

Full text of Hazard Statements: see section 16

SECTION 4: FIRST AID MEASURES

4.1. Description of first aid measures

General Advice	If symptoms persist, call a physician.
Eye Contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
Skin Contact	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
Ingestion	Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
Inhalation	Remove to fresh air. Get medical attention immediately if symptoms occur.
Self-Protection of the First Aider	No special precautions required.

4.2. Most important symptoms and effects, both acute and delayed

None reasonably foreseeable.

4.3. Indication of any immediate medical attention and special treatment needed

Notes to Physician	Treat symptomatically.
---------------------------	------------------------

SECTION 5: FIREFIGHTING MEASURES

5.1. Extinguishing media

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment. Water spray, carbon dioxide (CO₂), dry chemical, alcohol-resistant foam.

Extinguishing media which must not be used for safety reasons

No information available.

5.2. Special hazards arising from the substance or mixture

Thermal decomposition can lead to release of irritating gases and vapors.

Hazardous Combustion Products

None under normal use conditions.

5.3. Advice for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1. Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Use personal protective equipment as required.

6.2. Environmental precautions

Should not be released into the environment.

6.3. Methods and material for containment and cleaning up

Sweep up and shovel into suitable containers for disposal.

6.4. Reference to other sections

Refer to protective measures listed in Sections 8 and 13.

SECTION 7: HANDLING AND STORAGE

7.1. Precautions for safe handling

Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid ingestion and inhalation. Avoid contact with skin, eyes or clothing.

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Wash hands before breaks and after work.

7.2. Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Technical Rules for Hazardous Substances (TRGS) 510 Class 12
Storage Class (LGK) (Germany)

7.3. Specific end use(s)

Use in laboratories

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1. Control parameters

Exposure limits

List source(s): **EU** - Commission Directive (EU) 2019/1831 of 24 October 2019 establishing a fifth list of indicative occupational exposure limit values pursuant to Council Directive 98/24/EC and amending Commission Directive 2000/39/EC **UK** - EH40/2005 Work Exposure Limits, Fourth edition. Published 2020. **IRE** - 2021 Code of Practice for the Chemical Agents Regulations, Schedule 1. Published by the Health and Safety Authority

Component	The United Kingdom	European Union	Ireland
Ethyl alcohol	TWA: 1000 ppm TWA; 1920		STEL: 1000 ppm 15 min

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

	mg/m³ TWA WEL - STEL: 3000 ppm STEL; 5760 mg/m³ STEL		
Methyl alcohol	WEL - TWA: 200 ppm TWA; 266 mg/m³ TWA WEL - STEL: 250 ppm STEL; 333 mg/m³ STEL	TWA: 200 ppm 8 hr TWA: 260 mg/m³ 8 hr Skin	TWA: 200 ppm 8 hr. TWA: 260 mg/m³ 8 hr. STEL: 600 ppm 15 min STEL: 780 mg/m³ 15 min Skin

Biological limit values

List source(s):

Derived No Effect Level (DNEL) / Derived Minimum Effect Level (DMEL)

See table for values

Component	Acute effects local (Dermal)	Acute effects systemic (Dermal)	Chronic effects local (Dermal)	Chronic effects systemic (Dermal)
Ethyl alcohol 64-17-5 (<1.2)				DNEL = 343mg/kg bw/day
Dimethyl sulfoxide 67-68-5 (<0.6)				DNEL = 200mg/kg bw/day
Methyl alcohol 67-56-1 (<0.2)		DNEL = 20mg/kg bw/day		DNEL = 20mg/kg bw/day

Component	Acute effects local (Inhalation)	Acute effects systemic (Inhalation)	Chronic effects local (Inhalation)	Chronic effects systemic (Inhalation)
Ethyl alcohol 64-17-5 (<1.2)	DNEL = 1900mg/m³			DNEL = 950mg/m³
Dimethyl sulfoxide 67-68-5 (<0.6)			DNEL = 265mg/m³	DNEL = 484mg/m³
Methyl alcohol 67-56-1 (<0.2)	DNEL = 130mg/m³	DNEL = 130mg/m³	DNEL = 130mg/m³	DNEL = 130mg/m³
2-Methyl-2H-isothiazol-3-one 2682-20-4 (<0.01)	DNEL = 0.043mg/m³		DNEL = 0.021mg/m³	

Predicted No Effect Concentration (PNEC)

See values below.

Component	Fresh water	Fresh water sediment	Water Intermittent	Microorganisms in sewage treatment	Soil (Agriculture)
Ethyl alcohol 64-17-5 (<1.2)	PNEC = 0.96mg/L	PNEC = 3.6mg/kg sediment dw	PNEC = 2.75mg/L	PNEC = 580mg/L	PNEC = 0.63mg/kg soil dw
Dimethyl sulfoxide 67-68-5 (<0.6)	PNEC = 17mg/L	PNEC = 13.4mg/kg sediment dw		PNEC = 11mg/L	PNEC = 3.02mg/kg soil dw
Methyl alcohol 67-56-1 (<0.2)	PNEC = 20.8mg/L	PNEC = 77mg/kg sediment dw	PNEC = 1540mg/L	PNEC = 100mg/L	PNEC = 100mg/kg soil dw
2-Methyl-2H-isothiazol-3-one 2682-20-4 (<0.01)	PNEC = 3.39µg/L		PNEC = 3.39µg/L	PNEC = 0.23mg/L	PNEC = 0.0471mg/kg soil dw

Component	Marine water	Marine water sediment	Marine water intermittent	Food chain	Air
Ethyl alcohol 64-17-5 (<1.2)	PNEC = 0.79mg/L	PNEC = 2.9mg/kg sediment dw		PNEC = 0.38g/kg food PNEC = 0.72g/kg food	

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

Dimethyl sulfoxide 67-68-5 (<0.6)	PNEC = 1.7mg/L			PNEC = 0.7g/kg food	
Methyl alcohol 67-56-1 (<0.2)	PNEC = 2.08mg/L	PNEC = 7.7mg/kg sediment dw			
2-Methyl-2H-isothiazol-3-one 2682-20-4 (<0.01)	PNEC = 3.39µg/L		PNEC = 3.39µg/L		

8.2. Exposure controls

Engineering Measures

None under normal use conditions.

Personal protective equipment

Eye Protection

Wear safety glasses with side shields (or goggles) (European standard - EN 166)

Hand Protection

Protective gloves

Glove material	Breakthrough time	Glove thickness	EU standard	Glove comments
Disposable gloves	See manufacturers recommendations	-	EN 374	(minimum requirement)

Skin and body protection

Long sleeved clothing.

Inspect gloves before use.

Please observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves.
(Refer to manufacturer/supplier for information)

Ensure gloves are suitable for the task: Chemical compatability, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion.

Remove gloves with care avoiding skin contamination.

Respiratory Protection

No protective equipment is needed under normal use conditions.

Large scale/emergency use

Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced

Recommended Filter type: Particle filter

Small scale/Laboratory use

Maintain adequate ventilation

Environmental exposure controls

No information available.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1. Information on basic physical and chemical properties

Physical State

Liquid

Appearance

Odor

No information available

Odor Threshold

No data available

Melting Point/Range

No data available

Softening Point

No data available

Boiling Point/Range

No information available

Flammability (liquid)

No data available

Flammability (solid,gas)

Not applicable

Liquid

Explosion Limits

No data available

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

Flash Point	No information available	Method - No information available
Autoignition Temperature	No data available	
Decomposition Temperature	No data available	
pH	No information available	
Viscosity	No data available	
Water Solubility	No information available	
Solubility in other solvents	No information available	
Partition Coefficient (n-octanol/water)		
Component	log Pow	
Ethyl alcohol	-0.32	
Dimethyl sulfoxide	-1.35	
Methyl alcohol	-0.74	
Fluconazole	1.47	
2-Methyl-2H-isothiazol-3-one	-0.26	
Voriconazole	1.72	
Vapor Pressure	No data available	
Density / Specific Gravity	No data available	
Bulk Density	Not applicable	Liquid
Vapor Density	No data available	(Air = 1.0)
Particle characteristics	Not applicable (liquid)	

9.2. Other information

SECTION 10: STABILITY AND REACTIVITY

10.1. Reactivity

None known, based on information available

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous Polymerization

No information available.

Hazardous Reactions

None under normal processing.

10.4. Conditions to avoid

Incompatible products. Excess heat.

10.5. Incompatible materials

None known.

10.6. Hazardous decomposition products

None under normal use conditions.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Product Information

(a) acute toxicity;

Oral

Based on available data, the classification criteria are not met

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

Dermal
Inhalation

Based on available data, the classification criteria are not met
Based on available data, the classification criteria are not met

Toxicology data for the components

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Ethyl alcohol	LD50 = 7060 mg/kg (Rat)	-	20000 ppm/10H (Rat)
Dimethyl sulfoxide	LD50 = 28300 mg/kg (Rat)	LD50 = 40000 mg/kg (Rat)	LC50 > 5.33 mg/L (Rat) 4 h
Methyl alcohol	LD50 = 1187 – 2769 mg/kg (Rat)	LD50 = 17100 mg/kg (Rabbit)	LC50 = 128.2 mg/L (Rat) 4 h
Fluconazole	1271 mg/kg (Rat)	-	-
2-Methyl-2H-isothiazol-3-one	LD50 = 120 mg/kg (Rat) LD50 232 - 249 mg/kg (Rat)	LD50 = 200 mg/kg (Rabbit)	LC50 = 0.11 mg/L (Rat) 4 h
Voriconazole	100-300 mg/kg (Rat)	> 2000 mg/kg (Rat)	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-yl)methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-2-(1-methylpropyl)-	LD50 > 320 mg/kg (Rat)	-	-

(b) skin corrosion/irritation; No data available

(c) serious eye damage/irritation; No data available

(d) respiratory or skin sensitization;

Respiratory No data available
Skin No data available

Component	Test method	Test species	Study result
Methyl alcohol 67-56-1 (<0.2)	OECD Test Guideline 406 Guinea Pig Maximisation Test (GPMT)	guinea pig	non-sensitising

(e) germ cell mutagenicity; No data available

(f) carcinogenicity; No data available

The table below indicates whether each agency has listed any ingredient as a carcinogen

(g) reproductive toxicity; No data available

Component	Test method	Test species / Duration	Study result
Methyl alcohol 67-56-1 (<0.2)	OECD Test Guideline 416	Rat / Inhalation 2 Generation	NOAEC = 1.3 mg/l (air)

(h) STOT-single exposure; No data available

(i) STOT-repeated exposure; No data available

Target Organs None known.

(j) aspiration hazard; No data available

Symptoms / effects, both acute and No information available.

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

delayed

11.2. Information on other hazards

Endocrine Disrupting Properties Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

SECTION 12: ECOLOGICAL INFORMATION

12.1. Toxicity

Ecotoxicity effects

Component	Freshwater Fish	Water Flea	Freshwater Algae
Ethyl alcohol	Fathead minnow (Pimephales promelas) LC50 = 14200 mg/L/96h	EC50 = 9268 mg/L/48h EC50 = 10800 mg/L/24h	EC50 (72h) = 275 mg/l (Chlorella vulgaris)
Dimethyl sulfoxide	40 g/L LC50 96 h 33-37 g/L LC50 96 h	EC50 24h 7000 mg/L	EC50 96h 12350 - 25500 mg/L
Methyl alcohol	Pimephales promelas: LC50 > 10000 mg/L 96h	EC50 > 10000 mg/L 24h	
2-Methyl-2H-isothiazol-3-one	0.07 mg/l - 96h	0.18 mg/l - 48h	

Component	Microtox	M-Factor
Ethyl alcohol	Photobacterium phosphoreum: EC50 = 34634 mg/L/30 min Photobacterium phosphoreum: EC50 = 35470 mg/L/5 min	
Dimethyl sulfoxide	= 16000 mg/L EC50 Pseudomonas putida 16 h = 32 g/L EC50 Tetrahymena pyriformis 24 h = 77 mg/L EC50 Photobacterium phosphoreum 5 min	
Methyl alcohol	EC50 = 39000 mg/L 25 min EC50 = 40000 mg/L 15 min EC50 = 43000 mg/L 5 min	
2-Methyl-2H-isothiazol-3-one		10 (acute) 1 (chronic)

12.2. Persistence and degradability

No information available

Component	Degradability
Methyl alcohol 67-56-1 (<0.2)	DT50 ~ 17.2d >94% after 20d

12.3. Bioaccumulative potential

No information available

Component	log Pow	Bioconcentration factor (BCF)
Ethyl alcohol	-0.32	No data available
Dimethyl sulfoxide	-1.35	No data available
Methyl alcohol	-0.74	<10 dimensionless
Fluconazole	1.47	No data available
2-Methyl-2H-isothiazol-3-one	-0.26	No data available
Voriconazole	1.72	No data available

12.4. Mobility in soil

No information available

12.5. Results of PBT and vPvB assessment

No data available for assessment.

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

12.6. Endocrine disrupting properties

Endocrine Disruptor Information

12.7. Other adverse effects

Persistent Organic Pollutant
Ozone Depletion Potential

This product does not contain any known or suspected substance

This product does not contain any known or suspected substance

SECTION 13: DISPOSAL CONSIDERATIONS

13.1. Waste treatment methods

Waste from Residues/Unused Products

Waste is classified as hazardous. Dispose of in accordance with the European Directives on waste and hazardous waste. Dispose of in accordance with local regulations.

Contaminated Packaging

Dispose of this container to hazardous or special waste collection point.

European Waste Catalogue (EWC)

According to the European Waste Catalog, Waste Codes are not product specific, but application specific.

Other Information

Waste codes should be assigned by the user based on the application for which the product was used. Do not empty into drains.

SECTION 14: TRANSPORT INFORMATION

IMDG/IMO

Not regulated

14.1. UN number

14.2. UN proper shipping name

14.3. Transport hazard class(es)

14.4. Packing group

ADR

Not regulated

14.1. UN number

14.2. UN proper shipping name

14.3. Transport hazard class(es)

14.4. Packing group

IATA

Not regulated

14.1. UN number

14.2. UN proper shipping name

14.3. Transport hazard class(es)

14.4. Packing group

14.5. Environmental hazards

No hazards identified

14.6. Special precautions for user

No special precautions required

14.7. Maritime transport in bulk according to IMO instruments

Not applicable, packaged goods

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

SECTION 15: REGULATORY INFORMATION

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

International Inventories

X = listed, U.S.A. (TSCA), Canada (DSL/NDL), Europe (EINECS/ELINCS/NLP), Australia (AICS), Korea (KECL), China (IECSC), Japan (ENCS), Philippines (PICCS), Taiwan (TCSI), Japan (ISHL), New Zealand (NZIoC), Japan (ISHL). US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

Component	CAS No	EINECS	ELINCS	NLP	IECSC	TCSI	KECL	ENCS	ISHL
Human Serum	NA	-	-	-	-	-	-	-	-
Ethyl alcohol	64-17-5	200-578-6	-	-	X	X	KE-13217	X	X
Dimethyl sulfoxide	67-68-5	200-664-3	-	-	X	X	KE-32367	X	X
Methyl alcohol	67-56-1	200-659-6	-	-	X	X	KE-23193	X	X
Fluconazole	86386-73-4	-	-	-	-	X	-	-	X
2-Methyl-2H-isothiazol-3-one	2682-20-4	220-239-6	-	-	X	X	KE-24316	X	X
Voriconazole	137234-62-9	-	-	-	-	X	-	-	X
Posaconazole	171228-49-2	-	-	-	-	-	-	-	X
Isavuconazole	241479-67-4	-	-	-	-	-	-	-	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[4-[[2-(2,4-dichlorophenyl)- 2-(1H-1,2,4-triazol-1-yl)methyl]-1,3- dioxolan-4-yl]methoxy]phenyl]-1-pi peraziny]phenyl]-2,4-dihydro-2-(2- hydroxy-1-methylpropyl)-	112559-91-8	-	-	-	-	-	-	-	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[4-[[2-(2,4-dichlorophenyl)- 2-(1H-1,2,4-triazol-1-yl)meth yl]-1,3-dioxolan-4-yl]methoxy]phen yl]-1-piperaziny]phenyl]-2,4- dihydro-2-(1-methylpropyl)-	84625-61-6	-	-	-	-	X	-	-	X

Component	CAS No	TSCA	TSCA Inventory notification - Active-Inactive	DSL	NDL	AICS	NZIoC	PICCS
Human Serum	NA	-	-	-	-	-	-	-
Ethyl alcohol	64-17-5	X	ACTIVE	X	-	X	X	X
Dimethyl sulfoxide	67-68-5	X	ACTIVE	X	-	X	X	X
Methyl alcohol	67-56-1	X	ACTIVE	X	-	X	X	X
Fluconazole	86386-73-4	-	-	-	-	-	-	-
2-Methyl-2H-isothiazol-3-one	2682-20-4	X	ACTIVE	X	-	X	X	X
Voriconazole	137234-62-9	-	-	-	-	-	-	-
Posaconazole	171228-49-2	-	-	-	-	-	-	-
Isavuconazole	241479-67-4	-	-	-	-	-	-	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[4-[[2-(2,4-dichlorophenyl)- 2-(1H-1,2,4-triazol-1-yl)methyl]-1,3- dioxolan-4-yl]methoxy]phenyl]-1-pi peraziny]phenyl]-2,4-dihydro-2-(2- hydroxy-1-methylpropyl)-	112559-91-8	-	-	-	-	-	X	-
3H-1,2,4-Triazol-3-one, 4-[4-[4-[4-[[2-(2,4-dichlorophenyl)- 2-(1H-1,2,4-triazol-1-yl)meth yl]-1,3-dioxolan-4-yl]methoxy]phen yl]-1-piperaziny]phenyl]-2,4- dihydro-2-(1-methylpropyl)-	84625-61-6	-	-	-	-	-	-	-

Legend: X - Listed '-' - Not Listed

KECL - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

Authorisation/Restrictions according to EU REACH

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

Component	REACH (1907/2006) - Annex XIV - Substances Subject to Authorization	REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances	REACH Regulation (EC 1907/2006) article 59 - Candidate List of Substances of Very High Concern (SVHC)
Dimethyl sulfoxide	-	Use restricted. See item 75. (see link for restriction details)	-
Methyl alcohol	-	Use restricted. See item 69. (see link for restriction details)	-
2-Methyl-2H-isothiazol-3-one	-	Use restricted. See item 75. (see link for restriction details)	-

After the sunset date the use of this substance requires either an authorization or can only be used for exempted uses, e.g. use in scientific research and development which includes routine analytics or use as intermediate.

<https://echa.europa.eu/authorisation-list>

<https://echa.europa.eu/substances-restricted-under-reach>

<https://echa.europa.eu/candidate-list-table>

Component	CAS No	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements
Human Serum	NA	Not applicable	Not applicable
Ethyl alcohol	64-17-5	Not applicable	Not applicable
Dimethyl sulfoxide	67-68-5	Not applicable	Not applicable
Methyl alcohol	67-56-1	500 tonne	5000 tonne
Fluconazole	86386-73-4	Not applicable	Not applicable
2-Methyl-2H-isothiazol-3-one	2682-20-4	Not applicable	Not applicable
Voriconazole	137234-62-9	Not applicable	Not applicable
Posaconazole	171228-49-2	Not applicable	Not applicable
Isavuconazole	241479-67-4	Not applicable	Not applicable
3H-1,2,4-Triazol-3-one, 4-[4-[4-[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-2-(2-hydroxy-1-methylpropyl)-	112559-91-8	Not applicable	Not applicable
3H-1,2,4-Triazol-3-one, 4-[4-[4-[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-2-(1-methylpropyl)-	84625-61-6	Not applicable	Not applicable

Regulation (EC) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of dangerous chemicals

Not applicable

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work .

Take note of Directive 2000/39/EC establishing a first list of indicative occupational exposure limit values

Take note of Directive 94/33/EC on the protection of young people at work

Take note of Dir 92/85/EC on the protection of pregnant and breastfeeding women at work

National Regulations

UK - Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment

WGK Classification

Water endangering class = 1 (self classification)

Component	Germany - Water Classification (AwSV)	Germany - TA-Luft Class
Ethyl alcohol	WGK1	
Dimethyl sulfoxide	WGK1	
Methyl alcohol	WGK 2	

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

2-Methyl-2H-isothiazol-3-one	WGK3	
------------------------------	------	--

Component	France - INRS (Tables of occupational diseases)
Ethyl alcohol	Tableaux des maladies professionnelles (TMP) - RG 84
Dimethyl sulfoxide	Tableaux des maladies professionnelles (TMP) - RG 84
Methyl alcohol	Tableaux des maladies professionnelles (TMP) - RG 84

Component	Switzerland - Ordinance on the Reduction of Risk from handling of hazardous substances preparation (SR 814.81)	Switzerland - Ordinance on Incentive Taxes on Volatile Organic Compounds (OVOC)	Switzerland - Ordinance of the Rotterdam Convention on the Prior Informed Consent Procedure
Ethyl alcohol 64-17-5 (<1.2)		Group I	
Methyl alcohol 67-56-1 (<0.2)	Prohibited and Restricted Substances	Group I	

15.2. Chemical safety assessment

Chemical Safety Assessment/Reports (CSA/CSR) are not required for mixtures

SECTION 16: OTHER INFORMATION

Full text of H-Statements referred to under sections 2 and 3

H225 - Highly flammable liquid and vapor
H301 - Toxic if swallowed
H302 - Harmful if swallowed
H311 - Toxic in contact with skin
H314 - Causes severe skin burns and eye damage
H315 - Causes skin irritation
H317 - May cause an allergic skin reaction
H319 - Causes serious eye irritation
H331 - Toxic if inhaled
H335 - May cause respiratory irritation
H351 - Suspected of causing cancer
H360 - May damage fertility or the unborn child
H360D - May damage the unborn child
H361 - Suspected of damaging fertility or the unborn child
H370 - Causes damage to organs
H400 - Very toxic to aquatic life
H412 - Harmful to aquatic life with long lasting effects

Legend

CAS - Chemical Abstracts Service

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

IECSC - Chinese Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDL - Canadian Domestic Substances List/Non-Domestic Substances List

ENCS - Japanese Existing and New Chemical Substances

AICS - Australian Inventory of Chemical Substances

NZIoC - New Zealand Inventory of Chemicals

WEL - Workplace Exposure Limit

ACGIH - American Conference of Governmental Industrial Hygienists

DNEL - Derived No Effect Level

RPE - Respiratory Protective Equipment

TWA - Time Weighted Average

IARC - International Agency for Research on Cancer
Predicted No Effect Concentration (PNEC)

LD50 - Lethal Dose 50%

SAFETY DATA SHEET

Cascadion™ SM Antifungals Calibrator Set, Cascadion™ SM Antifungals Control 1,
Cascadion™ SM Antifungals Control 2, Cascadion™ SM Antifungals Control 3

Revision Date 07-Mar-2022

LC50 - Lethal Concentration 50%
NOEC - No Observed Effect Concentration
PBT - Persistent, Bioaccumulative, Toxic

EC50 - Effective Concentration 50%
POW - Partition coefficient Octanol:Water
vPvB - very Persistent, very Bioaccumulative

ADR - European Agreement Concerning the International Carriage of Dangerous Goods by Road
IMO/IMDG - International Maritime Organization/International Maritime Dangerous Goods Code
OECD - Organisation for Economic Co-operation and Development
BCF - Bioconcentration factor

ICAO/IATA - International Civil Aviation Organization/International Air Transport Association
MARPOL - International Convention for the Prevention of Pollution from Ships
ATE - Acute Toxicity Estimate
VOC - (Volatile Organic Compound)

Key literature references and sources for data

<https://echa.europa.eu/information-on-chemicals>
Suppliers safety data sheet, Chemadviser - LOLI, Merck index, RTECS

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

Physical hazards	On basis of test data
Health Hazards	Calculation method
Environmental hazards	Calculation method

Training Advice

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

Creation Date	07-Mar-2022
Revision Date	07-Mar-2022
Revision Summary	SDS sections updated.

**This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.
COMMISSION REGULATION (EU) 2020/878 amending Annex II to Regulation (EC) No
1907/2006**

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

End of Safety Data Sheet