

Occupational Safety and Health (Classification, Labelling and Safety Data Sheet of Hazardous Chemicals) Regulations 2013

**SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**
**Product Identifier**

Perihalan Produk:

Product Description:

Cat No. :

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

FE/1300/15; FE/1300/17; FE/1300/99V

**Relevant identified uses of the substance or mixture and uses advised against**

Recommended Use

Laboratory chemicals.

Uses advised against

No Information available

**Company**

 Thermo Fisher Scientific Fisher Scientific (M) Sdn Bhd  
 Hap Seng Business Park, Lot 01-03, 01-04 Aras 1 Unity Square,  
 No 12, Persiaran Perusahaan, Seksyen 23, 40300 Shah Alam,  
 Selangor Darul Ehsan, Malaysia.  
 Main line: +60 3-5525 7888

**Supplier**

E-mail address

Enquiry.my@thermofisher.com

Emergency Telephone Number

Tel: +03-5525 7888

CHEMTREC Malaysia 1-800-815-308 (Malay)

CHEMTREC Malaysia (Kuala Lumpur) +(60)-327884561 (Malay)

**SECTION 2: HAZARDS IDENTIFICATION**
**Classification of the substance or mixture**

|                                                    |                   |
|----------------------------------------------------|-------------------|
| Flammable liquids                                  | Category 2 (H225) |
| Acute oral toxicity                                | Category 3 (H301) |
| Acute dermal toxicity                              | Category 3 (H311) |
| Acute Inhalation Toxicity - Vapors                 | Category 3 (H331) |
| Specific target organ toxicity - (single exposure) | Category 1 (H370) |

**Label Elements**


Signal Word

Danger

Hazard Statements

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

H225 - Highly flammable liquid and vapor  
H370 - Causes damage to organs  
H301 + H311 + H331 - Toxic if swallowed, in contact with skin or if inhaled

## Precautionary Statements

### Storage

P405 - Store locked up

### Disposal

P501 - Dispose of contents/ container to an approved waste disposal plant

## Other Hazards

Toxic to terrestrial vertebrates  
This product does not contain any known or suspected endocrine disruptors

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

| Component      | CAS No  | Weight % |
|----------------|---------|----------|
| Methyl alcohol | 67-56-1 | 99.9     |
| Formic acid    | 64-18-6 | 0.1      |

## SECTION 4: FIRST AID MEASURES

### Description of first aid measures

#### Eye Contact

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.

#### Skin Contact

Wash off immediately with plenty of water for at least 15 minutes. Immediate medical attention is required.

#### Ingestion

Do NOT induce vomiting. Call a physician or poison control center immediately.

#### Inhalation

Remove to fresh air. Do not use mouth-to-mouth method if victim ingested or inhaled the substance; give artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper respiratory medical device. Immediate medical attention is required. If not breathing, give artificial respiration.

#### Self-Protection of the First Aider

Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination.

### Most important symptoms and effects, both acute and delayed

Difficulty in breathing. Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting.

### Indication of any immediate medical attention and special treatment needed

#### Notes to Physician

Treat symptomatically. Symptoms may be delayed.

## SECTION 5: FIREFIGHTING MEASURES

### Extinguishing media

Suitable Extinguishing Media

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

CO<sub>2</sub>, dry chemical, dry sand, alcohol-resistant foam. Water mist may be used to cool closed containers. Water mist may be used to cool closed containers.

## Extinguishing media which must not be used for safety reasons

No information available.

## Special hazards arising from the substance or mixture

Flammable. Risk of ignition. Vapors may form explosive mixtures with air. Vapors may travel to source of ignition and flash back. Containers may explode when heated. Vapors may form explosive mixtures with air.

## Hazardous Combustion Products

Carbon monoxide (CO), Formaldehyde.

## Advice for fire-fighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear. Thermal decomposition can lead to release of irritating gases and vapors.

## SECTION 6: ACCIDENTAL RELEASE MEASURES

### Personal Precautions, Protective Equipment and Emergency Procedures

Use personal protective equipment as required. Remove all sources of ignition. Evacuate personnel to safe areas. Keep people away from and upwind of spill/leak. Take precautionary measures against static discharges. Do not get in eyes, on skin, or on clothing.

### Environmental precautions

Should not be released into the environment.

### Methods and Material for Containment and Cleaning Up

Remove all sources of ignition. Soak up with inert absorbent material. Take precautionary measures against static discharges. Keep in suitable, closed containers for disposal. Use spark-proof tools and explosion-proof equipment.

### Reference to Other Sections

Refer to protective measures listed in Sections 8 and 13.

## SECTION 7: HANDLING AND STORAGE

### Precautions for Safe Handling

Use only under a chemical fume hood. Use spark-proof tools and explosion-proof equipment. Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharges. Do not breathe mist/vapors/spray. Do not get in eyes, on skin, or on clothing. Use only non-sparking tools. To avoid ignition of vapors by static electricity discharge, all metal parts of the equipment must be grounded.

### Conditions for Safe Storage, Including any Incompatibilities

Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from open flames, hot surfaces and sources of ignition. Flammables area. Keep away from heat, sparks and flame.

### Specific End Uses

Use in laboratories.

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

## Control Parameters

| Component      | Malaysia | ACGIH TLV                             | OSHA PEL                                                                                                                                                                                 |
|----------------|----------|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Methyl alcohol |          | TWA: 200 ppm<br>STEL: 250 ppm<br>Skin | (Vacated) TWA: 200 ppm<br>(Vacated) TWA: 260 mg/m <sup>3</sup><br>(Vacated) STEL: 250 ppm<br>(Vacated) STEL: 325 mg/m <sup>3</sup><br>Skin<br>TWA: 200 ppm<br>TWA: 260 mg/m <sup>3</sup> |
| Formic acid    |          | TWA: 5 ppm<br>STEL: 10 ppm            | (Vacated) TWA: 5 ppm<br>(Vacated) TWA: 9 mg/m <sup>3</sup><br>TWA: 5 ppm<br>TWA: 9 mg/m <sup>3</sup>                                                                                     |

| Component      | European Union                                               | The United Kingdom                                                                                               | Germany                                                                                                                                                                                                                                                  |
|----------------|--------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Methyl alcohol | TWA: 200 ppm 8 hr<br>TWA: 260 mg/m <sup>3</sup> 8 hr<br>Skin | WEL - TWA: 200 ppm TWA; 266 mg/m <sup>3</sup> TWA<br>WEL - STEL: 250 ppm STEL; 333 mg/m <sup>3</sup> STEL        | 100 ppm TWA MAK; 130 mg/m <sup>3</sup> TWA MAK Skin absorber                                                                                                                                                                                             |
| Formic acid    | TWA: 5 ppm (8hr)<br>TWA: 9 mg/m <sup>3</sup> (8hr)           | STEL: 15 ppm 15 min<br>STEL: 28.8 mg/m <sup>3</sup> 15 min<br>TWA: 5 ppm 8 hr<br>TWA: 9.6 mg/m <sup>3</sup> 8 hr | TWA: 5 ppm (8 Stunden). AGW - exposure factor 2<br>TWA: 9.5 mg/m <sup>3</sup> (8 Stunden). AGW - exposure factor 2<br>TWA: 5 ppm (8 Stunden). MAK<br>TWA: 9.5 mg/m <sup>3</sup> (8 Stunden). MAK<br>Höhepunkt: 10 ppm<br>Höhepunkt: 19 mg/m <sup>3</sup> |

## Exposure Controls

### Engineering Measures

Use only under a chemical fume hood. Use explosion-proof electrical/ventilating/lighting equipment. Ensure that eyewash stations and safety showers are close to the workstation location. Ensure adequate ventilation, especially in confined areas.

Wherever possible, engineering control measures such as the isolation or enclosure of the process, the introduction of process or equipment changes to minimise release or contact, and the use of properly designed ventilation systems, should be adopted to control hazardous materials at source

### Personal protective equipment

#### Eye Protection

Wear safety glasses with side shields (or goggles)

#### Hand Protection

Protective gloves

#### Skin and body protection

Long sleeved clothing

Inspect gloves before use.

Please observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information)

Ensure gloves are suitable for the task: Chemical compatibility, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion.

Remove gloves with care avoiding skin contamination.

#### Respiratory Protection

When workers are facing concentrations above the exposure limit they must use appropriate certified respirators

#### Recommended Filter type:

low boiling organic solvent Type AX Brown conforming to EN371

To protect the wearer, respiratory protective equipment must be the correct fit and be used and maintained properly

When RPE is used a face piece Fit Test should be conducted

## Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice

## Environmental exposure controls

No information available

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

|                                         |                           |                                             |
|-----------------------------------------|---------------------------|---------------------------------------------|
| Appearance                              | Clear, colorless solution |                                             |
| Physical State                          | Liquid                    |                                             |
| Odor                                    | Alcohol-like              |                                             |
| Odor Threshold                          | No data available         |                                             |
| pH                                      | 5                         |                                             |
| Melting Point/Range                     | No data available         |                                             |
| Softening Point                         | No data available         |                                             |
| Boiling Point/Range                     | 65 °C / 149 °F            | @ 760 mmHg                                  |
| Flash Point                             | 12 °C / 53.6 °F           | Method - No information available           |
| Evaporation Rate                        | No data available         |                                             |
| Flammability (solid,gas)                | Not applicable            | Liquid                                      |
| Explosion Limits                        | No data available         |                                             |
| Vapor Pressure                          | No data available         |                                             |
| Vapor Density                           | No data available         | (Air = 1.0)                                 |
| Specific Gravity / Density              | 0.79                      |                                             |
| Bulk Density                            | Not applicable            | Liquid                                      |
| Water Solubility                        | No information available  |                                             |
| Solubility in other solvents            | No information available  |                                             |
| Partition Coefficient (n-octanol/water) |                           |                                             |
| Component                               | log Pow                   |                                             |
| Methyl alcohol                          | -0.74                     |                                             |
| Formic acid                             | -1.9                      |                                             |
| Autoignition Temperature                | 455 °C / 851 °F           |                                             |
| Decomposition Temperature               | No data available         |                                             |
| Viscosity                               | No data available         |                                             |
| Explosive Properties                    |                           | Vapors may form explosive mixtures with air |
| Oxidizing Properties                    | No information available  |                                             |

## SECTION 10: STABILITY AND REACTIVITY

### Reactivity

None known, based on information available.

### Chemical Stability

Stable under normal conditions.

### Possibility of Hazardous Reactions

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

## Hazardous Polymerization Hazardous Reactions

Hazardous polymerization does not occur.  
None under normal processing.

## Conditions to Avoid

Incompatible products. Heat, flames and sparks. Keep away from open flames, hot surfaces and sources of ignition.

## Incompatible Materials

Strong oxidizing agents. Strong acids. Acid anhydrides. Acid chlorides. Strong bases. Metals. Peroxides.

## Hazardous Decomposition Products

Carbon monoxide (CO). Formaldehyde.

## SECTION 11: TOXICOLOGICAL INFORMATION

### Information on Toxicological Effects

#### Product Information

##### (a) acute toxicity;

|            |            |
|------------|------------|
| Oral       | Category 3 |
| Dermal     | Category 3 |
| Inhalation | Category 3 |

### Toxicology data for the components

| Component      | LD50 Oral                      | LD50 Dermal                   | LC50 Inhalation               |
|----------------|--------------------------------|-------------------------------|-------------------------------|
| Methyl alcohol | LD50 = 1187 – 2769 mg/kg (Rat) | LD50 = 17100 mg/kg ( Rabbit ) | LC50 = 128.2 mg/L ( Rat ) 4 h |
| Formic acid    | LD50 = 1100 mg/kg ( Rat )      | -                             | LC50 = 7.85 mg/L ( Rat ) 4 h  |

(b) skin corrosion/irritation; Based on available data, the classification criteria are not met

(c) serious eye damage/irritation; Based on available data, the classification criteria are not met

##### (d) respiratory or skin sensitization;

|             |                                                                  |
|-------------|------------------------------------------------------------------|
| Respiratory | Based on available data, the classification criteria are not met |
| Skin        | Based on available data, the classification criteria are not met |

| Component                          | Test method                                                       | Test species | Study result    |
|------------------------------------|-------------------------------------------------------------------|--------------|-----------------|
| Methyl alcohol<br>67-56-1 ( 99.9 ) | OECD Test Guideline 406<br>Guinea Pig Maximisation Test<br>(GPMT) | guinea pig   | non-sensitising |

(e) germ cell mutagenicity; Based on available data, the classification criteria are not met  
Mutagenic effects have occurred in experimental animals

(f) carcinogenicity; Based on available data, the classification criteria are not met  
There are no known carcinogenic chemicals in this product

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

**(g) reproductive toxicity;** Based on available data, the classification criteria are not met

| Component                          | Test method             | Test species / Duration          | Study result              |
|------------------------------------|-------------------------|----------------------------------|---------------------------|
| Methyl alcohol<br>67-56-1 ( 99.9 ) | OECD Test Guideline 416 | Rat / Inhalation<br>2 Generation | NOAEC =<br>1.3 mg/l (air) |

**(h) STOT-single exposure;** Category 1

**Results / Target organs** Optic nerve, Central nervous system (CNS).

**(i) STOT-repeated exposure;** Based on available data, the classification criteria are not met

**Target Organs** None known.

**(j) aspiration hazard;** Based on available data, the classification criteria are not met

**Other Adverse Effects** The toxicological properties have not been fully investigated.

**Symptoms / effects, both acute and delayed** Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting.

**Endocrine Disrupting Properties** Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

## SECTION 12: ECOLOGICAL INFORMATION

### Ecotoxicity effects

| Component      | Freshwater Fish                               | Water Flea            | Freshwater Algae   | Microtox                                                                        |
|----------------|-----------------------------------------------|-----------------------|--------------------|---------------------------------------------------------------------------------|
| Methyl alcohol | Pimephales promelas:<br>LC50 > 10000 mg/L 96h | EC50 > 10000 mg/L 24h |                    | EC50 = 39000 mg/L 25 min<br>EC50 = 40000 mg/L 15 min<br>EC50 = 43000 mg/L 5 min |
| Formic acid    | Leuciscus idus: LC50 =<br>46-100 mg/L/96h     | EC50 = 34 mg/L/48h    | EC50 = 25 mg/L/96h | EC50 = 46.7 mg/L/17h                                                            |

**Persistence and degradability** No information available

**Persistence** Persistence is unlikely, based on information available.

| Component                          | Degradability                  |
|------------------------------------|--------------------------------|
| Methyl alcohol<br>67-56-1 ( 99.9 ) | DT50 ~ 17.2d<br>>94% after 20d |

**Bioaccumulative potential** Bioaccumulation is unlikely

| Component      | log Pow | Bioconcentration factor (BCF) |
|----------------|---------|-------------------------------|
| Methyl alcohol | -0.74   | <10 dimensionless             |
| Formic acid    | -1.9    | 0.22 dimensionless            |

**Mobility in soil** The product contains volatile organic compounds (VOC) which will evaporate easily from all surfaces. Will likely be mobile in the environment due to its volatility. Disperses rapidly in air.

**Endocrine Disruptor Information** This product does not contain any known or suspected endocrine disruptors

# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

Other adverse effects No information available

## SECTION 13: DISPOSAL CONSIDERATIONS

### Waste treatment methods

**Waste from Residues/Unused Products**

Waste is classified as hazardous Dispose of in accordance with the European Directives on waste and hazardous waste Dispose of in accordance with local regulations

### **Contaminated Packaging**

Dispose of this container to hazardous or special waste collection point. Empty containers retain product residue, (liquid and/or vapor), and can be dangerous Keep product and empty container away from heat and sources of ignition

### **Other Information**

Waste codes should be assigned by the user based on the application for which the product was used Do not flush to sewer Can be landfilled or incinerated, when in compliance with local regulations

## SECTION 14: TRANSPORT INFORMATION

### IMDG/IMO

UN-No UN1230  
Hazard Class 3  
Subsidiary Hazard Class 6.1  
Packing Group II  
Proper Shipping Name METHANOL SOLUTION

### Road and Rail Transport

UN-No UN1230  
Hazard Class 3  
Subsidiary Hazard Class 6.1  
Packing Group II  
Proper Shipping Name METHANOL SOLUTION

### IATA

UN-No UN1230  
Hazard Class 3  
Subsidiary Hazard Class 6.1  
Packing Group II  
Proper Shipping Name METHANOL SOLUTION

**Special Precautions for User** No special precautions required

## SECTION 15: REGULATORY INFORMATION

### Safety, health and environmental regulations/legislation specific for the substance or mixture

**International Inventories** X = listed

| Component      | EINECS    | TSCA | DSL | PICCS | ENCS | ISHL | IECSC | AICS | KECL     |
|----------------|-----------|------|-----|-------|------|------|-------|------|----------|
| Methyl alcohol | 200-659-6 | X    | X   | X     | X    | X    | X     | X    | KE-23193 |
| Formic acid    | 200-579-1 | X    | X   | X     | X    | X    | X     | X    | KE-17233 |



# SAFETY DATA SHEET

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

| Component      | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification | Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements | Rotterdam Convention (PIC) | Basel Convention (Hazardous Waste) |
|----------------|-------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|----------------------------|------------------------------------|
| Methyl alcohol | 500 tonne                                                                                 | 5000 tonne                                                                               |                            |                                    |
| Formic acid    |                                                                                           |                                                                                          |                            | Annex I - Y34                      |

## National Regulations

**Persistent Organic Pollutant**  
**Ozone Depletion Potential**

This product does not contain any known or suspected substance  
This product does not contain any known or suspected substance

## SECTION 16: OTHER INFORMATION

### Legend

**CAS** - Chemical Abstracts Service

**EINECS/ELINCS** - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**IECSC** - Chinese Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**ENCS** - Japanese Existing and New Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**NZIoC** - New Zealand Inventory of Chemicals

**WEL** - Workplace Exposure Limit

**ACGIH** - American Conference of Governmental Industrial Hygienists

**RPE** - Respiratory Protective Equipment

**LC50** - Lethal Concentration 50%

**POW** - Partition coefficient Octanol:Water

**TWA** - Time Weighted Average

**IARC** - International Agency for Research on Cancer

**LD50** - Lethal Dose 50%

**EC50** - Effective Concentration 50%

**ADR** - European Agreement Concerning the International Carriage of Dangerous Goods by Road

**IMO/IMDG** - International Maritime Organization/International Maritime Dangerous Goods Code

**OECD** - Organisation for Economic Co-operation and Development

**BCF** - Bioconcentration factor

**ICAO/IATA** - International Civil Aviation Organization/International Air Transport Association

**MARPOL** - International Convention for the Prevention of Pollution from Ships

**ATE** - Acute Toxicity Estimate

**VOC** - (Volatile Organic Compound)

### Key literature references and sources for data

<https://echa.europa.eu/information-on-chemicals>

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

**Revision Date**

23-Mar-2025

**Revision Summary**

Not applicable.

**In accordance with local and national regulations: Occupational Safety and Health (Classification, Labelling and Safety Data Sheet of Hazardous Chemicals) Regulations 2013**

### Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

# **SAFETY DATA SHEET**

Methanol with 0,1% (v/v) Formic acid for HPLC Gradient analysis

Revision Date 23-Mar-2025

---

**End of Safety Data Sheet**